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INTEGRATED DILATOMETER-TENSIOMETER TO MEASURE THE STATE OF TENSION OF GLAZED BODIES

par / by Chiara Venturelli, Mariano Paganelli, Expert System Solutions, Modène, Italie / Modena, Italy

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Some frequent glazing defects such as crazing are caused by an inadequate state of tension between the glaze and the ceramic body, which is strongly related to various parameters. The most important are the coefficients of thermal expansion of glaze and body, their coupling temperature and the softening temperature of the glaze. In this work, a novel optical technique is proposed with the purpose to investigate the state of tension of glazed bodies and prevent the compatibility defects mentioned above.

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The authors want to thank the Department of Engineering of Materials and Environment of the University of Modena and Reggio Emilia for the kind collaboration.

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There are several methods for measuring the dimensional changes of the sample with a beam of light and they can be summarised in two categories: reflected beam and direct beam.

Reflected beam

The first optical dilatometer was invented by Abbe and Fizeau in the second half of the XVIIIth century. This design uses a reflected beam of monochromatic light and the measurement of the displacement is carried out by counting the interference fringes between the forward going beam and the reflected beam. But in order to measure thermal expansion up to high temperatures the sample has to be set on a sample holder inside the furnace and the sample holder also has its own thermal expansion. In order to achieve good accuracy it is necessary to measure the expansion of the sample holder and subtract it from the actual expansion of the specimen. The best way to do this is to split the laser beam into two beams of light, which are reflected by the top of the sample and by the sample holder or by both ends of the specimen. This method proved to be very accurate, but it is limited by the reflectance of the surface of the specimen. If the specimen is non-reflective, or becomes nonreflective during the test, it is then necessary to use a mirror. Working in this way, this method loses the advantage of being non-contact and becomes substantially similar to the electronic dilatometer.

Direct beam

The direct beam measuring system overcomes this problem because the specimen is measured by the image that it projects on a CCD when it is lighted by a direct beam of light. Using a beam of light with a short wavelength (for example a blue light with a wavelength below 0.5 µm) and a very high resolution CCD (actual resolution of 0.5 µm per pixel of the camera) it is possible to achieve good resolution.

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During firing, a glaze applied on a tile undergoes some transformations, which can be simply listed as: evaporation of the forming water of the clay components, glass transition and softening. In correspondence of this point, the glazed part of the glaze starts to melt, giving

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as a result a continuous liquid layer. During cooling, the glaze viscosity increases until the glaze becomes rigid and it starts to contract simultaneously with the ceramic body. Since a glass is much more resistant to compression than to traction, the thermal contraction of the glaze during cooling should not exceed the thermal contraction of the ceramic body. If the glaze is in a state of traction with respect to the ceramic body, the defect of crazing may occur and fractures in the vitreous layer appear. On the other hand, peeling is the defect, which may take place if the compressive state is too high. This defect is however less frequent if single firing is adopted, but also in this case some other problems can be present. For example, in the tile industry, a common defect is the fact that the tiles may result difficult to cut with a traditional scoring tile cutter, or the glaze may chip off during the scoring. Another problem for large size ceramic tiles is that an excessive state of compression of the glaze leads to convex tiles. The presence of these defects can be however not simply impute to the coefficients of thermal expansion mismatch between the glaze and the ceramic body, because another parameter should be considered: the coupling temperature between the glaze and the ceramic body. It is important to underline, in fact, that below the coupling temperature, glaze and body are bonded together and they contract in a different amount because of the difference in their thermo-mechanical behaviour, not simply because of the difference in their thermal expansion coefficients. Also glazes with coefficients of thermal expansion lower than the body may give rise to problems of crazing, if the coupling temperature is not right.

OLD METHOD

The old method used to measure the state of tension of a glaze consisting of

an analysis performed with two different instruments:


- a traditional dilatometer, which allows to obtain the thermal expansion curves of the glaze and the body and their coefficient of thermal expansion (C.T.E). The common laboratory practice is to carry out thermal expansion tests up to 400-500 °C. This is due to the fact that traditional electronic dilatometers work in contact with the sample to analyse: if the material softens or develops liquid phases, as in the case of glazes, the pressure imposed by the mechanical system deforms the material and the measuring system risks to be damaged;
- the Steger tensiometer, which determines the coupling temperature (in general, it is much higher than 400-500 °C, limit of temperature reached with the dilatometric tests). The test performed with the old Steger tensiometer (see fig. 1) is quite complex and presents the disadvantage that the specimen to be tested has to be prepared on purpose in laboratory and cannot be cut directly from an actual tile. The sample must be very long (30 cm) and it is not placed completely inside the kiln because its extremities must be external in order to do the measurement. As a consequence, the sample temperature is not homogeneous.

A consideration that should be done is that it is useless to perform dilatometric tests up to a temperature lower than the glaze-body coupling point. In fact, in this case the data obtained from the thermal expansion tests cannot be combined together with the Steger tensiometer tests results, in order to determine the state of tension of a glazed body. A technique able to investigate the thermal behaviour of materials over their softening point is required.

A NOVEL OPTICAL TECHNIQUE

The technique of analysis presented in this work is based on an apparatus (Misura Flex-ODLT, produced by Expert System Solutions and shown in figure 3),

Figure 7 – Courbes de dilatation thermique de l'émail et du carreau
Thermal expansion curves of glaze and body

which combines the functions of the two instruments mentioned above. This instrument works with no contact with the samples tested: thus, their behaviour is never modified by the measuring system. The specimens can be cut directly from an industrial tile; then they are placed into the kiln, which is equipped with an automatic temperature controller and can reach a heating rate of 30 °C per minute and a maximum temperature of 1,400 °C. All the experimental part of this work was performed by using the optical instrument shown in fig. 3.

EXPERIMENTAL PART

The first experimental step was the determination of the coupling temperature between glaze and body. The specimen is a bar of glazed material which has been cut directly from a tile industrially fired. It is not necessary to prepare a specific specimen, made on purpose for the test, like in the old Steger tensiometer. Two holding rods, which are 70 mm spaced, support the sample (see fig. 4). A camera frames the centre of the sample, which moves downward or upward during the test.

In fig. 5 the flexion of the glazed body is shown as a function of the temperature. It is possible to see that the sample is bending downward up to 600 °C and then it starts to bend upward. The negative flexion is due to the fact that the body has a higher coefficient of thermal expansion with respect to the C.T.E of the glaze. The flexion reaches the maximum in correspondence to the quartz transition, when the body shows the maximum coefficient of thermal expansion. After this point, the flexion starts to decrease, because the reduction of the viscosity of the glassy phase allows the


Figure 8 – Les courbes de dilatation thermique de l'émail et du carreau recouvrent la courbe de flexion du carreau émaillé (émail en compression) Thermal expansion curves of glaze and body overlapped with the flexion curve of the glazed body (glaze in compression)

the whole cooling phase below the coupling temperature. For this reason, the glaze remains shorter with respect to the body when cooled down. The simple comparison of the two coefficients of thermal expansion is misleading: nothing anomalous can be seen.

In fact, the glaze coefficient of thermal expansion ($5,38 \cdot 10^{-6} \text{ }^\circ\text{C}^{-1}$) is lower with respect to the body one ($6,92 \cdot 10^{-6} \text{ }^\circ\text{C}^{-1}$).

The problem lays in the coupling temperature, which is found in correspondence of the sharp change in slope in the glaze thermal expansion. Finally, the eventual presence of an engobe layer is a further variable which has to be taken into consideration and can be studied with this experimental procedure.

CONCLUSIONS

Many advantages can be found by using the new optical technique presented in this work. First of all, the specimens to be tested can be obtained directly from an industrial glazed product: this was not possible with the old Steger tensiometer. The state of tension of glazed bodies can be determined by using only the new instrument Misura, which combines both the dilatometer and the tensiometer.

Finally, it is possible to carry out thermal expansion and flexion tests even on materials that cannot be measured with traditional dilatometers, because of the softening or the developing of a liquid phase. For this reason, the thermal expansion test results can be correlated with the Steger tensiometer test results, in order to prevent the most serious compatibility defect, i.e. crazing.

viscous flow and absorbs completely the tensions. At the beginning of the cooling phase, the glaze is liquid and follows the body contraction without developing tensions. Then, an abrupt variation of inclination in the flexion curve occurs: this means the glaze has become rigid again, starting to build up tension and causing the flexion of the sample. The coupling temperature is identified in correspondence of this point: in our case is 700 °C.

The second experimental phase, performed with the same instrument, consists of the thermal expansion tests performed on the fired glaze and ceramic body. In this case, two beams of light illuminate both the ends of a specimen 50 mm long placed horizontally into the furnace and two digital cameras capture the images of the last few hundreds microns of each tip. The specimen, completely free to expand or contract, is measured by the image that it projects on a CCD.

In fig. 7, the thermal expansion curves of the glaze and the ceramic body are represented. In the case of the glaze thermal expansion curve, a comparison with the results provided by a traditional electronic push rod dilatometer could be done. The first part of the curve will be very similar, but more significant differences will appear at higher temperature (over the softening point). The curve obtained with the optical instrument Misura shows a shortening of the sample, but in this case the sample is actually becoming shorter because the surface tension is rounding the edges, not because the push rod is deforming it. In fig. 8, the flexion and the two thermal expansion curves are overlapped. The

glaze thermal expansion curve is translated so that it coincides with the body thermal expansion curve in correspondence of the coupling temperature. In this way, the two curves do not coincide anymore at the origin (room temperature). This difference is indeed the level of traction/compression, which has established between the glaze and the ceramic body. In the case studied, it is possible to see that the ceramic body after cooling is shorter than the glaze and therefore keeps it in compression. The tension conditions of the glazed body are far away to lead to problems of crazing.

The opposite situation is provided in fig. 9, obtained following the method used for the previous example. In this case, the glaze may present problems of crazing, because it is in a state of traction during


Figure 9 – Les courbes de dilatation thermique de l'émail et du carreau recouvrent la courbe de flexion du carreau émaillé (émail en traction) Thermal expansion curves of glaze and body overlapped with the flexion curve of the glazed body (glaze in traction)