Optical non contact dilatometry

by Mariano Paganelli

Expert System, Modena (Italy)

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The study of the behaviour and properties of ceramic materials during the heating process is one of the most critical steps in the development of innovative products in the fields of advanced and traditional ceramic materials.

With the optical non-contact dilatometer, the range of research can be extended to temperatures higher than the dilatometric softening point because the measuring system does not alter the sample's behaviour while performing the measurement.

This makes it possible to faithfully reproduce the industrial thermal cycle and to monitor the sample's behaviour even during the sintering or expansion phases through a precise measurement of the size variation of the material under investigation.

The thermomechanical behaviour of ceramic materials is usually studied up to temperatures slightly higher than the dilatometric softening temperature through the use of electronic dilatometers with a linear displacement transducer (LVDT).

These transducers can easily reach resolution levels higher than 1/10000 of a millimetre but they work using a mechanical system that keeps the core of the sensor in contact with the sample.

Mechanical contact is achieved using a small bar made of refractory material which is kept at a slight pressure on the sample by a spring or a leverage.

The slight pressure exerted by the measuring system does not disturb the measurement of the thermomechanical properties as long as the sample remains rigid, but when the sample begins to soften the measurement is no longer valid since the slight pressure of the measurement system modifies the sample's behaviour.

The non-contact optical dilatometer solves this problem completely since the behaviour of the sample is not modified by the measuring system.

In spite of this major advantage, however, we must

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remember that the resolution of an optical system is limited by the wavelength of the light used to perform the measurement.

Since visible light has a wavelength spectrum that varies from 0.3 to 0.7 microns, in optimal conditions using blue light a resolution of 3/10000 of a millimetre can be achieved.

For commercial instruments operating in filtered visible light we can consider the standard resolution as 5/10000 of a millimetre, that is five times lower than an electronic instrument using displacement sensors such as the LVDT.

This kind of dilatometer has numerous practical applications in all fields where it is important to ascertain the thermomechanical behaviour of the material even at temperatures where the sample displays a viscous behaviour

The field of traditional and advanced ceramic materials is perhaps the ideal application for this kind of measuring instrument, since all these materials have a wide range of temperatures over which their mechanical behaviour ranges from purely elastic to viscous-elastic and in the end becomes primarily viscous.

For example, glass has an essentially elastic behaviour at temperatures considerably lower than the glass transition temperature (T_g), a visco-elastic behaviour in the temperature interval around T_g , and finally a primarily viscous behaviour at temperatures considerably higher than T_g .

Ceramic bodies behave in a primarily rigid manner up to temperatures close to the sintering point and then assume a viscous-elastic behaviour with a viscosity that decreases exponentially with increasing temperature.

Listed below are some typical examples of applications on ceramic materials: a traditional glaze frit, a glass-ceramic material, a body for porous single-firing and a non-conventional body that forms cristobalite during firing.

A ceramic frit

The thermal expansion curve obtained using the optical dilatometer on a sample of ceramic glaze frit is very similar to that obtained with an electronic dilatometer (Fig. 1).

The measured values are similar for the interval up to the glass transition temperature (T_g) (segment A-B).

The curve segment beyond T_g (where the glass has an expansion coefficient characteristic of the

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liquid phase since it is in an undercooled liquid state), has the same slope but extends to higher temperatures than those measured using an electronic dilatometer (segment B-C₁-C).

The dilatometric softening temperature is considerably higher and the descending curve segment beyond the dilatometric softening temperature (segment C-D) has a negative slope with a much lower value than that measured with a mechanical instrument (segment C₁-D₁).

The differences observed using a non-contact instrument are all due to the fact that the measuring system does not exert pressure on the sample. In practical terms, information measured at low temperatures is basically identical although the optical instrument makes the measurements with a lower resolution.

A glass-ceramic material

Glass-ceramic materials are characterised by the fact that they are obtained in glassy phase and then transformed into crystal phases through a suitable heat treatment.

The study of the transition between the glass and crystalline phase is far from simple and requires specific equipment. The instruments typically used for this kind of research are X-ray diffraction (XRD) and differential thermal analysis (DTA). The electronic dilatometer is not conventionally used to study the transition between the glass and crystalline phases since the transition occurs at temperatures higher than T_g where the viscosity is sufficiently low to permit the

growth of the crystal nuclei. This temperature level is usually higher than the dilatometric softening temperature and is therefore beyond the scope of traditional measuring instruments.

The optical dilatometer on the other hand provides a new possibility for research in this field. Since it does not exert any pressure on the sample, it can monitor the process at higher temperatures than electronic dilatometers.

The example shown in figure 2 represents the behaviour of a glass-ceramic material sample.

The initial part of the curve, up to the dilatometric softening temperature, is entirely similar to that measured with an electronic dilatometer, but the second part of the curve shows the real behaviour of material

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heated at 30°C/min. After a descending segment in which the glass undergoes an apparent contraction due to a reduction in viscosity and consequent softening of the sample, a phase of marked expansion takes place. Crystallisation phenomena within the vitreous mass ($T_{\text{crystallisation}} = 910^{\circ}\text{C}$) make the material rigid again and it passes from the typically viscous behaviour of the glass phase at temperatures higher than the glass transition temperature to the elastic behaviour of the newly formed crystalline phases.

This kind of information may be invaluable in the development of new materials, since the actual expansion of the material in the crystallisation phase can give rise to deformation phenomena which otherwise could not be understood.

A body for porous single-firing

Bodies for porous single-firing are characterised by an ability to reach a good level of mechanical strength in spite of the fact that their density after firing is lower than their density prior to firing. This is because at a certain temperature solid state reactions begin and trigger the growth of newly-formed crystalline phases which are in turn responsible for variations in the material's mechanical properties.

The example shown in figure 3 shows the behaviour of a body for porous single-firing subjected to a temperature gradient of 60°C/min., which corresponds to the heating rate used in the industry for a firing cycle of 30 minutes.

The curve shows that the body undergoes a shrinkage phase between 900°C and 1000°C (segment B-C), then expansion up to 1150°C (segment C-D). The segment D-C represents cooling of the sample, showing that after firing the samples return to a size just slightly smaller than the initial size of the dry sample.

This curve was obtained with a dilatometric resolution of 6/10000 of a millimetre over the entire temperature interval.

The most useful piece of information we can derive from an analysis of this graph is certainly the growth interval of the newly-formed crystalline phases, which extends over the entire segment C-D of the curve. This kind of body has its highest crystalline growth at 1125°C (point D) when subjected to this heating gradient

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
Fig. 3. Sintering curve of a porous single-firing body recorded during firing and cooling, obtained with dilatometric resolution (resolution 0.6 microns, sample length 15 mm) using the optical dilatometer Misura 3

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(60°C/min.). At this temperature there is a small degree of shrinkage compared to the initial size of the dry sample. This is the temperature at which the highest level of crystallinity is obtained, hence the greatest stability against moisture-induced expansion.

A special ceramic material

Figure 4 shows the behaviour of a special ceramic material with a very high expansion coefficient. The upper portion of the curve shows the behaviour of a raw sample subjected to a heating cycle with 60°C/min. gradient, followed by the behaviour of the fired sample during natural cooling down to ambient temperature. In the first part of the graph we can observe a complex sintering process causing the material's crystalline structure to be entirely modified.

At 600°C a sintering reaction begins, which then continues up to 1020°C, followed by an expansion phase due to recrystallisation of the glass phases.

The expansion curve of the fired sample displays the presence of a large quantity of cristobalite (transition at 200°C) which has formed during the firing cycle and was not present in the initial material.

With this kind of test, the sample's progress can be monitored during the firing cycle without altering its behaviour, after which the test can be continued on the same sample during the cooling phase.

Conclusions

The availability of a laboratory instrument that can make thermomechanical measurements without contact with the sample and with high heating rates opens up new areas of research which until now have been beyond the scope of research and development laboratories working in the traditional ceramic sector.

A very large amount of information can certainly be obtained from non-contact measurements, especially in the field of innovative materials whose behaviour is complex and difficult to understand.

This new opportunity to monitor the processes of sintering, recrystallisation and growth of neoformation phases on site and during the thermal cycle is a powerful tool in the hands of researchers working in the field of traditional and advanced ceramic materials.

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